



Auscultation Manikin Update Package

KERI™ — LF04024U

GERI™ — LF04004U

Instruction Manual



Life/form® Products by Nasco



Congratulations for choosing the **GERI™/KERI™** Auscultation Update Package. This simple design makes the manikin easy to use and care for. A diagnosis of the condition selected by the instructor can be made by comparing the variation in sounds occurring at different sites — just like a real person! Lung sounds can be detected at five anterior and 12 posterior locations, and students can practice auscultation at six anterior heart sites. Before using this manikin with your students, please familiarize yourself with this manual.

List of Components

- **GERI™/KERI™** Manikin Update
- SmartScope™
- Remote Control



Figure 1

General Instructions for Use

The **GERI™/KERI™** Update comes ready to be attached to your existing manikin. To do this, you will need to use the arms, legs, and head of your existing manikin. The arms can be

removed by bending the elbows backward approximately 90° to align the keyholes at the shoulders. (**See figure 1.**) Then pull outward, away from the shoulder of the manikin. Reverse this procedure to attach the arms to the new **GERI™/KERI™** Auscultation Torso.



Figure 2

The legs detach from the body at the hips by bending the legs backward approximately 150° so that the feet are near the shoulders and the keyholes are aligned. (**See figure 2.**) Then pull outward, away from the body. Reverse this procedure to attach to the new torso.



Figure 3

The head may be removed by rotating it 180° backward to align the keyholes, and then pulling up. The internal stomach reservoir will come out of the body through the neck opening. Reverse this procedure to attach the head to the new torso.

Gastrostomy Care

Lavage and Gavage

A flanged hole simulating an abdominal incision for the insertion of a feeding tube is included on the upper torso for performing lavage and gavage. Inside the upper torso is a reservoir (bag) with a maximum capacity of 500 cc, which is attached to the gastrostomy opening with a two-part coupler. The bag is removed by pushing in on the clear L-shaped button and pulling it straight away. To attach, simply push the coupler body (with bag attached) onto the coupler insert (part with a black O-ring) that is protruding from inside the torso. You will hear a slight snap when the connection is complete. With extensive use, these two coupler parts may no longer snap, so the clear L-shaped button will have to be manually pulled all the way out after connecting the two parts. Test the connection by pulling on the coupler body to ensure that it is locked. This gastrostomy feature is designed for use with a 16 French feeding tube. It is recommended that the tip of the feeding tube be lubricated before inserting. Make sure the reservoir is straight and flat before attempting to simulate actual feeding with a liquid.

About the Auscultation Simulator

The **GERI™/KERI™** Auscultation Manikin duplicates heart and lung conditions selected by the instructor by wireless remote control. The student should palpate to identify the correct auscultation sites. They will hear different heart and lung sounds as the SmartScope™ is moved from site to site.

The simulator has six heart sites and seven lung sound sites on the anterior surface, with 10 locations on the posterior surface and midaxillary sites. The remote control can select from 12 different heart conditions, as well as 16 lung conditions. The instructor can select any condition easily, enabling the student to compare sounds and make a diagnosis.

The remote control does not have to be pointed directly at the manikin or stethoscope to operate. One remote control will operate multiple sets of SmartScopes™ and manikins simultaneously. Great for group instruction. The range of the remote control is up to 100 feet.

General Instructions for Use

To begin using the Auscultation Trainer, press the red power button. This turns on the remote control and sends a signal to activate the stethoscope as well. After the unit is activated, the display will be in the “status” mode, displaying the current menu settings for the heart and lung conditions.

To select a new condition, press either the heart or lung button. This will put the display into the menu mode. The user can select a condition by either using the number buttons, or by viewing the conditions in sequence using the scroll button. After the condition is selected, press the enter button to activate.

The remote control can be programmed to these heart conditions:

- 01 Normal
- 02 Aortic regurgitation
- 03 Pulmonary stenosis
- 04 Mitral stenosis
- 05 Holosystolic
- 06 Mid-systolic
- 07 S3 Gallop
- 08 S4 Gallop
- 09 Systolic click
- 10 Atrial septal defect
- 11 PDA
- 12 VSD

The remote control can be programmed to these lung conditions:

- 01 Normal Lung
- 02 Normal Vesicular
- 03 Wheezes
- 04 Mono wheeze
- 05 Fine crackle
- 06 Coarse crackle
- 07 Ronchi crackle
- 08 Stridor
- 09 Cavernous
- 10 Bronchovesicular
- 11 Bronchial
- 12 Pulmonary edema
- 13 Infant
- 14 Friction rub
- 15 Egophony
- 16 Pectoriloquy

Consult the enclosed laminated card for lung and heart sound site locations.

To listen to the selected sounds, place the earpieces in ears, angled in a forward position. Place the diaphragm of the SmartScope™ over one of the appropriate sites on the manikin (see diagram on back).

As an alternative, an Amplifier/Speaker System (LF01189U) can be used, allowing an entire classroom to hear at the same time. To connect the speaker, plug the speaker cord into the speaker jack on the top of the Smart Scope™ box. When a speaker is connected to the SmartScope™, the earpieces will not work. Place the diaphragm of the SmartScope™ over the appropriate sites on the manikin.

Note: The remote has a battery saver that shuts the unit down after eight minutes if the remote is left on the same setting. To prevent it from shutting down, select a different heart or lung sound within the eight-minute period of time.

Available Supplies


LF01144U Additional SmartScope™

LF01148U Replacement Remote Control

LF01189U Amplifier/Speaker System

GERI™/KERI™ THE NURSING SKILLS MANIKIN

Component Replacement Parts



Other Available *Life/form*® Simulators

LF00698U	Adult Injectable Arm (Light)	LF01174U	NG Tube & Trach Skills
LF00855U	Male Catheterization	LF01184U	Venatech IM & Sub Q
LF00856U	Female Catheterization	LF01193U	Special Needs Baby
LF00901U	Prostate Examination	LF03000U	CPARLENE ® Series
LF00906U	Ostomy Care	LF03601U	Adult Airway Management Trainer with Stand
LF00929U	Surgical Bandaging	LF03602U	Adult Airway Management Manikin
LF00957U	Enema Administration	LF03609U	Child Airway Management Trainer with Stand
LF00958U	Pediatric Injectable Arm	LF03616U	Child CRiSis ™ Manikin
LF00961U	Intramuscular Injection	LF03617U	Deluxe Child CRiSis ™ Manikin with Arrhythmia Tutor
LF00984U	Breast Examination	LF03620U	PALS Update Kit
LF00995U	Arterial Puncture Arm	LF03623U	Infant Airway Management Trainer with Stand
LF00999U	Pediatric Injectable Head	LF03632U	Child Intraosseous Infusion/ Femoral Access Leg on a Stand
LF01005U	First Aid Arm	LF03633U	Child Airway Management Trainer Torso
LF01008U	Intradermal Injection Arm	LF03693U	Basic Buddy ® CPR Manikin
LF01012U	Heart Catheterization (TPN)	LF03699U	"Airway Larry" Airway Management Trainer
LF01019U	Ear Examination	LF03709U	Infant CRiSis ™ Manikin
LF01027U	Peritoneal Dialysis	LF03720U	Baby Buddy ™ Infant CPR Manikin
LF01028U	Suture Practice Arm	LF03750U	Bariatric CPR Manikin
LF01034U	Suture Practice Leg	LF03770U	Chest Tube
LF01036U	Spinal Injection	LF03953U	CRiSis ™ Manikin, Complete
LF01037U	Hemodialysis Practice Arm	LF03955U	Deluxe CRiSis ™ Manikin
LF01038U	Episiotomy Suturing Set	LF03956U	Deluxe "Plus" CRiSis ™ Manikin
LF01042U	Suture Kit	LF03965U	Adult CRiSis ™ Auscultation Manikin
LF01062U	Pelvic, Normal & Abnormal	LF03966U	Adult CRiSis ™ Auscultation Manikin with ECG Simulator
LF01063U	Stump Bandaging, Upper	LF04000U	GERI ™ / KERI ™ Manikin Series
LF01064U	Stump Bandaging, Lower	LF04200U	Adult Sternal Intraosseous Infusion
LF01069U	Cervical Effacement	LF06001U	CPR Prompt® Adult/Child Manikin
LF01070U	Birthing Station	LF06012U	CPR Prompt® Infant Manikin
LF01082U	Cricothyrotomy	LF06200U	CPR Prompt® Keychain Rescue Aid
LF01083U	Tracheostomy Care	LF06204U	CPR Prompt® Rescue and Practice Aid
LF01084U	Sigmoidoscopic Examination		
LF01087U	Central Venous Cannulation		
LF01095U	Blood Pressure Arm		
LF01108U	Infant Intraosseous Infusion		
LF01121U	Advanced IV Arm		
LF01131U	Venipuncture and Injection Arm		
LF01139U	Advanced IV Hand		
LF01142U	Auscultation Trainer		
LF01143U	Testicular Exam		
LF01152U	Male & Female Catheter		
LF01155U	Advanced CPR Dog		
LF01162U	Venatech IV Trainer		

Nasco Fort Atkinson

901 Janesville Avenue, P.O. Box 901
Fort Atkinson, Wisconsin 53538-0901

1.800.558.9595

eNasco.com • E-mail: lifeform@eNasco.com

COPYRIGHT © NASCO 2001 PRINTED IN U.S.A. NP 98-01/RV 10-15