
INSTRUCTION MANUAL

Economy Adult Sani-Manikin 100-2144





Assembly and Instructions

The parts needed to assemble the manikin are:

(1) manikin torso, (1) manikin head (1) airway sample, (2) pine tree push pins, (2) washers and (2) hairpin clips. (1) Rubberband

STEPS:

- 1. Insert the pin into the ear and through the neck, and then add the washer and connecting pin on the inside of the manikin's head.
- 2. Insert the pine tree push pins (2) into the back of the manikin. (See Figure 1)




Figure 1

- **3**. Place the rubber band over the pine tree push pins.
- 4. Fold the closed end of the airway lengthwise four times to simplify pushing the airway through the manikin's mouth. Tilt the head back and insert the closed end of the airway into and through the back of the tube inside the manikin's mouth.

(See Figure 2)




Figure 2

- **5.** Pull the closed end of the airway down the back of the manikin, leaving approximately 2 inches remaining outside of the manikin's mouth to use as a face shield.
- **6.** Flatten the closed end of the airway against the manikin's back and secure it under the rubber band.

(See Figure 3)



Figure 3

7. Spread the open end of the airway over the manikin's nose and chin. Finger sweep to open the bag to ensure there are no twists and that the airway is open. When the individual completes airway inflation, have them remove the bag so the next student can install their airway. The airway can be retained

Assembly and Instructions

by individual for use during training sesson, but it is necessary to discard the airway at the end of the session. (See Figure 4)

remove the han on the next etudent can inetall their




Figure 4

GENERAL INSTRUCTONS FOR USE

To ventilate the manikin use rescue breathing or a bag valve mask. An open airway and adequate ventilations will expand the airway, lifting the manikin up off the floor. To eliminate cross contamination the airway is designed for a single user. Airways are to be disposed of at the end of the training session. Reduce adhesive build up from hands free pads by applying a thin coat of automotive wax to chest locations.

CLEANING AND MAINTAINING ECONOMY SANI CPR MANIKIN

Remove used airway and dispose. Clean manikin with warm soapy water. Dry the manikin before storing. Do not immerse in fluids. Store in a dry location to avoid absorbing moisture.

WARRANTY

Nasco Healthcare warrants this product to be free from defects in materials and/or workmanship for a period of three years from the date of purchase, as evidenced by the date on the invoice of the product shipment to the end user. The warranty expressly does not cover abuse, accidental or purposeful damage, or any form of modification to the product. Nasco Healthcare reserves the right to either replace affected parts or the entire unit, at their sole discretion, after investigating and reviewing the actual product and the damage. In most instances, a digital photo of the product in question showing the damage will help qualify a product for return to the factory. At no time will any product be accepted at the factory without proper return authorization issued by Nasco Healthcare.



Nasco Healthcare 16 Simulaids Drive Saugerties, NY 12477 1-833-NASCOHC (627-2642) info@nascohealthcare.com www.nascohealthcare.com