



ALS Cardiac Training Model - Basic



- Simulator for adult resuscitation training and evaluation
- Lifelike body structure with silicon for realistic training
- Various emergency simulations available
- Rechargeable battery with completely wireless environment via Wi-Fi and Bluetooth communication
- Pupillary reflex; patient assessment and evaluation
- Generate related pulse and check ECG change from patient monitor during chest compressions
- Simultaneous chest compressions and ventilation
- Palpate pulse at eight sites (left radial and brachial, left and right carotid/femoral/dorsal)
- Change the pulse strength, and turn each pulse off/on
- Virtual patient monitor layout can be adjusted by selecting required items
- When running a prepared scenario, the current patient status can be jumped to the next status by using a transition
- Scenario and checklist can be exported or imported
- Result data can be printed out by using a printer connected to a laptop
- Detailed CPR results and log can be saved in Excel format
- The event log can be sorted by item or time order



01530 512425



shop@simulaids.co.uk




SimulaidsTM

ALS Cardiac Training Model - Basic



Storage case



Full-body manikin power adaptor



Laptop for AppFacilitator



Wi-Fi router



All-in-one PC for AppMonitor



Sphygmomanometer



Simulated blood




Lubricant



Bottle and tube for IO



Bottle and tube for IV



Syringe for cleaning IV/IO



IO module (tibia)



IO/IM module (humerus)




IV pad



Defibrillation electrode cover for lead wires



Defibrillation electrode cover for pads



⚡ Defibrillation

● ECG

● Pulse

● Heart sound

● Lung sound



01530 512425



shop@simulaids.co.uk



/simulaidsltd