

BT Inc BP/Pulse Assessment Simulator (BT-CEAB2)

- Realistic Practice of B/P examination with a stethoscope
- Setting the Systolic/Diastolic BP (BP range: 35~240Hg)
- Adjusting the pulse rate per minute (Pulse rate range: 40~140-times/min)
- Full control of Systolic/Diastolic B/P & pulse rate
- Selecting Auscultatory gap (On/Off)
- Choosing pulse intensity (3 stages)
- Adjusting Korotkoff sound's volume (5 stages)
- Pulse palpation at radial & brachial arteries.
- Auscultate Korotkoff's sound of brachial artery
- And more...



01530 512425



shop@simulaids.co.uk



Simulaids™



SimulaidsTM