



BT Inc Catheterisation Simulator (BT-CSCM / BT-CSCF)

- Life-like practical training
- Positioning similar to a real patient
- Any water-soluble lubricant for clinical application can be applied
- Training with various sizes of catheters is available
- Realistic practice with attachable extended legs
- Non-return valve: Prevents leaking after taking the catheter out
- Genital organs and buttocks are detachable for easy maintenance
- Foreskin (detachable) (M)
- Battery powered alarm sounds when a catheter inserted into the vagina (F)



01530 512425



shop@simulaids.co.uk



Simulaids™

