



CENTURIS PRO BED

CUSTOMER VALUE

The **Centuris Pro** bed is designed to provide simplicity, robustness and compliance to safety standards. The intuitive operation helps to support ease of adaptability for healthcare professionals in short stay / medical and surgical locations. As part of the **Hill-Rom 900** family with proven reliability and standardised parts, helps to make maintenance efficient. A solution that can meet the needs when under strict budget management.



KEEPING PATIENTS SAFE

- Intuitive egress grip for patient turning and getting out of bed
- Low height indicator
- Steam clean validation



IMPACTING OUTCOMES

- **SlideGuard** technology designed to help reduce compression, friction and shear
- Compatible with Baxter surfaces
- Head of bed angle indicator
- Ready to be used thanks to trusted reliability:
 - Material strength simulation
 - Design Validation strong process
 - Quality control
 - Ease of maintenance



OPTIMISING WORKFLOW

- **Auto Contour** and **SlideGuard** technologies designed to help minimise patient migration in bed
- Integrated X-Ray Accessory
- Patient and caregiver controls clearly separated

TECHNICAL SPECIFICATIONS

Safe working load		220 kg
Bed weight		137 kg
Low position		39.7 cm
SlideGuard technology sliding backrest		up to 12 cm
Protection class		IPX4
Sliderail height		39.3 cm
Bed angles	Head section	0 to 65°
	Thigh section	0 to 28°
	Foot section	-22° to -3°
	Trendelenburg/Reverse Trendelenburg	-17 to +17°
Sleep deck (HPL and metal slatted)	Width	90 cm
	Length	200 cm
Casters		150 mm Integral

STANDARD **CENTURIS PRO** BED PORTFOLIO FEATURES INCLUDE

- Split siderails
- Mobile head section
- 150 mm castors
- Linen holder
- Bed extension
- **Line-Of-Site** head of bed angle indicator
- Brake at foot end
- Steering wheel at headend
- Bilateral emergency CPR control
- Removable head and foot boards
- Removable HPL and metal slatted sleep deck
- Validated for steam cleaning

WHAT IS INCLUDED INTO THE **CENTURIS PRO**



Head of bed angle indicator



Trendelenburg &
Reverse Trendelenburg



Soft drop, ergonomically designed
side rails



Hand pendant control



Linen Holder



Integrated bed extension

The Centuris Pro bed is intended for use in acute, general and ambulatory care or care during long hospital stays for adult patients with a weight ≥ 40 kg, height ≥ 146 cm and BMI ≥ 17 . Medical Device (Regulation 2017/745): Class I

Manufacturer: Hill-Rom S.A.S. – Z.I. du Talhouet, 56330 Pluvigner, France.

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Baxter recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Baxter reserves the right to make changes without notice in design, specifications and models. The only warranty Baxter makes is the express written warranty extended on the sale or rental of its products.

baxterhealthcare.co.uk

Baxter Healthcare Ltd
Wallingford Road, Compton, Newbury, Berkshire, RG20 7QW

Baxter, Auto Contour, Centuris Pro, Hill-Rom 900, Line-of-Site and SlideGuard are trademarks of Baxter International Inc. or its subsidiaries.

UKI-CS136-250001 [v1.0] 02/2025