

Child Airway Management Trainer with Stand



SKU: LF03609

A simulated head of a 5-year-old child featuring realistic size and anatomy.



- ✓ Training stand to mount head
- ✓ Pump spray lubricant
- ✓ Hard carrying case

Perfect for practicing airway management skills on paediatric patients.

Realistic anatomy includes mouth, tongue, oral pharynx, larynx, epiglottis, vocal cords, trachea, and oesophagus.

Separate lungs for auscultation. Inflatable stomach bladder indicates oesophageal insufflation, manual carotid pulse. Practice suctioning and oral, nasal, and intubation techniques.



Simulaids