



Injectable Training Model

LF00952U & LF00953U

INSTRUCTION CARD



ABOUT THE SIMULATOR

Nasco developed a lifelike trainer to enhance training for intradermal, subcutaneous, and intramuscular injecting. The material feels like real human flesh. It can be pinched and punctured for injections repeated times.

GENERAL INSTRUCTIONS FOR USE

Dry injecting is recommended. The trainer will reject water or any liquid injected into the trainer. Solvents or corrosive materials will damage the simulator.

SN

Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.

Note: Save this card should you need to reference the Serial Number above.

5-Year Warranty

CARE OF THE SIMULATOR

Normal soil can be removed with mild soap and water. Nasco Cleaner will remove stubborn stains. Apply Nasco Cleaner to the soiled area and wipe clean with a soft cloth or paper towel. Allow to dry thoroughly before storing. Never place the simulator on newsprint or any kind of printed paper, plastic, or ballpoint pen. These materials will transfer indelible stains.

AVAILABLE SUPPLIES

LF01049U	Suture Kit with DVD Collection
LF07000U	Suture Replacement
LF07000(A)U	Suture Replacement, Pkg. of 5
LF09919U	Nasco Cleaner
LF03461U	Optional Carry Case
C00879U	Mayo Hegar Needle Holder

Other Available *Life/form* Simulators

LF00890U	Interactive Suture Trainer, White	LF01042U	Suture Kit
LF00892U	Interactive Suture Trainer, Black	LF01046U	Facial Suturing Module Set
LF00961U	Intramuscular Injection Simulator	LF01047U	Pediatric Suture Head Kit
LF01008U	Intradermal Injection Simulator	LF01131U	Venipuncture & Injection Demonstration Arm
LF01028U	Suture & Stapling Practice Arm	LF01184U	Venatech IM & SubQ Simulator, White
LF01034U	Suture & Stapling Practice Leg	LF01185U	Venatech IM & SubQ Simulator, Black
		WA18238U	Diabetic Injection Pad
		WA28395U	Injection Belly



ONLINE: NascoHealthcare.com
E-MAIL: lifeform@enasco.com
PHONE: 1.800.431.4310
ADDRESS: 901 Janesville Ave.
Fort Atkinson, WI 53538