

Simulaids® Infant Airway Management Trainer



SKU: 101-115

Designed to teach oral and nasal endotracheal intubation procedures on Infants.



- ✓ 1 x Completely Assembled Infant Airway Trainer
- ✓ Lubricant
- ✓ Carry bag
- ✓ Instructions

An excellent choice for student training. Comes with teeth, uvula, vocal cords, glottis, epiglottis, larynx, arytenoid cartilage, trachea, esophagus, and inflatable lungs and stomach.

Colour highlighted vocal cords for easy laryngoscope viewing. Accepts 2.5 mm endotracheal tubes.

Mounted on practice board for easy class use.



Simulaids