

TERi™ Androgynous Geriatric Patient Trainer

A comprehensive male/female elderly patient care trainer for medical simulation. Evolving from our GERi™/KERi™ full-body manikins, it has a more realistic look and feel, accurately representing the human anatomy. Weight is distributed to represent a real patient, for lifting and carrying.

Each version is affordable and durable with upgradable items.

- 1. Geriatric Patient Care Trainer** – General patient care, daily living assistance simulation platform.
- 2. Geriatric Patient Skills Trainer** – Physical skills practice simulation platform with injection/IV, catheterisation, cardiovascular and respiratory support.
- 3. Geriatric Patient Simulator** – Comprehensive and complete simulation platform including patient monitoring capability.

Nasogastric Tube Procedures^{2,3}

Placement, lavage, gavage, feeding, and suctioning.



Grooming and Daily Living Assistance^{1,2,3}

Bathing, clothing changes, denture placement/removal, hearing aid placement/removal, eye irrigation and optic drops, ostomy and urinary catheter care, and hair care.



PROM: Passive Range of Motion^{1,2,3}

Seating, laying, and finger/toe manipulation with no effort from patient.



Cardiovascular Support^{2,3}

Carotid pulses, CPR with quality performance metrics and auscultation.



Indwelling Male and Female Urinary Catheterisation^{2,3}

Includes support for urethral and suprapubic catheters.



Advanced Skills^{2,3}

Urinary catheterisation, enema administration, pelvic and prostate examination.



Intravenous and Subcutaneous Injection Sites^{2,3}, optional on 1

Multiple locations, equipped with standard arm IM injection pads.



Respiratory Support^{2,3}

Manual ventilation with chest rise, oxygen administration, tracheostomy care, and airway/trach suctioning. Airway intubation, ventilation and auscultation.



Enhanced Aesthetics^{1,2,3}

Soft silicone skin with noticeable normal and cancerous skin moles, stage 1 sacral ulcer and reddened skin folds.



This manikin is intended to be used as an instructional aid for teaching purposes. Actual product may vary slightly from photo.



01530 512425



shop@simulaidstld.co.uk



Simulaidstld™

Product Name	TERi™ Androgynous Geriatric Patient Trainer
Needs/Curricula	General capability carefully selected to support practice on: <ul style="list-style-type: none"> ▪ Essentials of patient care ▪ Movement assistance and fall prevention ▪ Grooming and daily living assistance ▪ Basic and advanced nursing skills ▪ Cardiovascular and respiratory emergencies Suitable for long-term facilities and emergency room training
Configuration	<ul style="list-style-type: none"> ▪ Light or medium skin tones ▪ Configurable to fit unique training needs and budgets ▪ A range of accessories available to enhance the baseline version, maximising versatility, training benefit, and return on investment
Accessories	IV Arm – LF04303 BP Arm – LF04304
Consumables	<ul style="list-style-type: none"> ▪ IV Arm <ul style="list-style-type: none"> - Skins - Veins kit - IM pads - Fluid supply bag ▪ IM pads for thighs and buttocks ▪ Urinary catheterisation: <ul style="list-style-type: none"> - Male genitalia foreskins - Foley catheter - Urinary catheterisation bladders ▪ Lung bags ▪ Stomach bags ▪ Stoma bags ▪ Gastrostomy bags ▪ Simulated blood ▪ Lubricant
Product Code/SKU	Geriatric Patient Care Trainer (Light) - LF04300 Geriatric Patient Skills Trainer (Light) - LF04301 Geriatric Patient Simulator (Light) - LF04302
Dimensions/Weight	Dimensions: Five feet Ship Weight: 65lbs



01530 512425



shop@simulaid.co.uk



/simulaidsttd